



Silicone Implants

NEW STERILE FACIAL IMPLANTS

Hanson Medical has recently designed and manufactured a new line of formerly patented facial implants featuring superb balance and symmetry. These implants, similar to some you may already be using, are offered at a 25% to 40% savings over our competitors. They are provided sterile with free sizer sets offered complimentary with your initial order of each style. For your custom implant needs, we offer free moulage modeling kits and will send you a model of your design for approval prior to manufacturing.

ABOUT US

Hanson Medical's principal staff members have been serving the Plastic and Reconstructive Surgery Industry for over 25 years. Over the years we have worked with many of the masters to design, manufacture, and market quality facial and body contouring silicone elastomer implants. Our High Quality Manufacturing Operation, located in Paso Robles, California, is certified by the International Standards Organization. We are compliant with the FDA Code of Federal Regulations for Good Manufacturing Practices. We have passed every inspection by the FDA and the California Food and Drug Branch in excellent fashion. Our implant grade raw materials are manufactured by Nusil Corporation of Carpinteria California. We offer world class implants and friendly efficient customer service.

TABLE OF CONTENTS

Solid Silicone Implants

Malar Implants	3-6
Chin Implants	7-12
Nasal Implants	12-13
Body Contouring Implants	14



HANSON
MEDICAL, INC.

BIOCOMPATIBILITY STUDIES

Hanson Medical has performed a vast battery of biocompatibility studies on our finished implants. These studies include Cytotoxicity, Mutagenicity, Pyrogenicity, Chromosomal DNA Aberration Testing, Muscle Implantation Studies, Teratogenicity and others.

SYMMETRIC IMPLANTS

Recent research has discovered that the main key to facial attractiveness is symmetry. The female models shown all have been gifted with facial symmetry. We have made a concentrated effort to create implants with near perfect symmetry in order to create a beautifully balanced surgical outcome for your patients. Our implants may also be scalpel adjusted by the surgeon for the patient lacking ideal facial symmetry.

CUSTOM IMPLANTS

Hanson Medical offers rapid custom implant fabrication with an average turn time of one week. We will send you a model of your implant upon request for your final adjustment and approval. Custom implant design kits are complimentary with your custom implant purchase.

STERILE IMPLANTS


Our entire line of silicone elastomer facial implants are provided sterile. We are the first company to have validated a hydrogen peroxide gas sterilization process. Our unique processes have eliminated concerns with gamma radiation and ethylene oxide.

APPRECIATION

We thank you for reviewing this brochure and for considering Hanson Medical as your new partner in Plastic Surgery.

Thank you for operating with us.

**HANSON
MEDICAL, INC.**




Ideal Symmetry



MALAR STYLE 1

Anatomical Design

CAT #	Size	A	B	C
SM1-1	S	3.9	1.95	0.35
SM1-2	M	4.4	2.25	0.40
SM1-3	L	5.1	2.7	0.50
M1S	Malar Style 1 Sizer Set			




For the subtle augmentation of the cheek area. Extends along the Zygoma laterally.



MALAR STYLE 2

Malar Shield

CAT #	Size	A	B	C
SM2-1	S	4.95	2.4	0.33
SM2-2	M	5.15	2.8	0.40
SM2-3	L	5.25	3.2	0.42
SM2-4	XL	5.65	3.6	0.50
M2S	Malar Style 2 Sizer Set			




Creates a rounder cheek contour with a greater surface area. Compare to *Implantech® TMS Series



MALAR STYLE 3

Oval Design

CAT #	Size	A	B	C
SM3-1	S	3.8	1.9	0.40
SM3-2	L	4.4	2.1	0.45
M3S	Malar Style 3 Sizer Set			



This implant is designed for very subtle augmentation of the cheek.



HANSON
MEDICAL, INC.

to order please call +91 91155 55555 or mail info@mitocon.in




MALAR STYLE 4

Extended Malar Shield

CAT #	Size	A	B	C
SM4-1	S	4.3	2.7	0.5
SM4-2	M	4.7	2.85	0.55
SM4-3	L	5.35	3.25	0.6
SM4-4	XL	5.55	3.6	0.7
M4S	Malar Style 4 Sizer Set			

Creates a rounded sweeping contour of the cheek. It extends over the infra malar area to provide downward augmentation.




MALAR STYLE 5

Anatomical High Projection

CAT #	Size	A	B	C
SM5-1	S	4.5	2.8	0.52
SM5-2	M	4.9	3.0	0.60
SM5-3	L	5.55	3.3	0.72
M5S	Malar Style 5 Sizer Set			

For maximum malar projection. It also augments the lateral orbital rim and the zygoma.




MALAR STYLE 6

Brennan Midfacial Filler

CAT #	Size	A	B	C
SM6-1	S	4.6	1.6	0.7
SM6-2	M	4.6	1.6	0.9
SM6-3	L	4.6	1.6	1.1
M6S	Malar Style 6 Sizer Set			

This implant is concave convex and sits securely in the infra malar sulcus. Designed by George Brennan MD



to order please call +91 91155 55555 or mail info@mitocon.in

HANSON
MEDICAL, INC.






MALAR STYLE 7

Inframalar

CAT #	Size	A	B	C
SM7-1	S	4.0	1.6	0.45
SM7-2	M	4.4	1.8	0.50
SM7-3	L	4.9	2.0	0.55
SM7-4	XL	5.3	2.2	0.62
M7S	Malar Style 7 Sizer Set			




Fills the inframalar sulcus and augments the mid cheek area. Compare to *Implantech® BMS Series.



MALAR STYLE 9

Tear Sulcus Anatomical Design

CAT #	Size	A	B	C
SM9-1	S	4.25	2.2	0.35
SM9-2	M	4.10	2.1	0.40
SM9-3	L	4.75	2.38	0.50
M9S	Malar Style 9 Sizer Set			




This implant is designed to fill in the tear groove below the eye. Brightens the cheeks in patients with flat infraorbital areas. Compare to *Implantech® FTT Series.



MALAR STYLE 10

Extended Tear Sulcus

CAT #	Size	A	B	C
SM10	S	5.8	2.6	0.51
M10S	Malar Style 10 Sizer			




Designed for wider augmentation of the infraorbital area. Compare to *Implantech® EFTT Series.



HANSON
MEDICAL, INC.


to order please call +91 91155 55555 or mail info@mitocon.in



MALAR STYLE 11

Inframalar II

CAT #	Size	A	B	C
SM11-1	S	4.0	1.9	0.43
SM11-2	M	4.4	2.22	0.46
SM11-3	L	4.9	2.30	0.45
SM11-4	XL	5.3	2.55	0.53
M11S	Malar Style 11 Sizer Set			




Provides projection for the midface above and below the infra malar sulcus. Compare to *Implantech® BMSII Series.



MALAR STYLE 12

Dual Projection

CAT #	Size	A	B	C
SM12-1	S	4.5	2.8	0.45
SM12-2	M	5.2	3.1	0.46
SM12-3	L	5.6	3.5	0.60
SM12-4	XL	5.95	3.7	0.63
M12S	Malar Style 12 Sizer Set			




Compare to *Implantech® CSM Series.



CARVING BLOCKS

Various Durometers

CAT #	A	B	C
BL13	7.5	15.5	1.3
BL20	7.5	15.5	2.0



Durometer choices: 5-Very Soft, 10-Soft, 15-Medium Soft, 20-Medium, 30-Firm
Specify color: Fleshstone, Clear or White

to order please call +91 91155 55555 or mail info@mitocon.in

HANSON
MEDICAL, INC.






CHIN STYLE 1

Anatomical

CAT #	Size	A	B	C	D
SCH1-1	S	5.6	1.1	0.6	8.7
SCH1-2	M	5.8	1.3	0.7	10.4
SCH1-3	L	6.0	1.4	0.8	10.8
SCH1-4	XL	6.1	1.4	0.9	10.8

CH1S Chin Style 1 Sizer Set

May be augmented with CH17



This chin provides forward and lateral projection. The smaller sizes are ideal for female patients.




CHIN STYLE 2

Flat Back Chin

CAT #	Size	A	B	C	D
SCH2-1	XS	5.7	1.1	0.6	8.5
SCH2-1.7	S	5.5	1.1	0.7	8.7
SCH2-2	M	6.6	1.3	1.0	9.6
SCH2-3	L	7.0	1.4	1.0	9.7
SCH2-4	XL	6.7	1.45	1.2	9.7

CH2S Chin Style 2 Sizer Set

Designed by Kurt Wagner, MD



The style 2 chin is flat backed and offers substantial projection in the larger sizes. The smaller sizes are perfect for female patients.




CHIN STYLE 3

Chin Jowl Style I

CAT #	Size	A	B	C	D	E
SCH3-1	S	5.7	1.2	0.6	0.55	10.8
SCH3-2	M	5.9	1.3	0.7	0.6	11.0
SCH3-3	L	6.0	1.35	0.9	0.5	11.0
SCH3-XL		5.8	1.4	1.05	0.7	11.0

CH3S Chin Style 3 Sizer Set



Chin Jowl Style I Offers increased projection to fill out the jowl groove.



HANSON
MEDICAL, INC.


to order please call +91 91155 55555 or mail info@mitocon.in



CHIN STYLE 4

Flat Back Button Chin

CAT #	Size	A	B	C
SCH4-1	S	3.2	1.1	0.5
SCH4-2	M	3.6	1.5	0.6
SCH4-3	L	3.8	1.5	0.6
SCH4-4	XL	4.3	1.8	0.7
CH4S	Chin Style 4 Sizer Set			




Designed for the patient desiring forward projection only. Ideal for use accompanying facelifts and necklifts.



CHIN STYLE 5

Concave Back Button Chin

CAT #	Size	A	B	C
SCH5-1	S	4.0	1.3	0.40
SCH5-2	M	4.0	1.3	0.60
SCH5-3	L	4.0	1.3	0.65
SCH5-4	XL	4.0	1.5	0.7
SCH5-5	XXL	4.7	1.75	0.8
CH5S	Chin Style 5 Sizer Set			




Offers finer feathered edges than the Chin Style 4 (SCH 4).



CHIN STYLE 6

KD Cleft Chin

CAT #	Size	A	B	C	D
SCH6-1	S	5.6	1.3	1.3	8.7
SCH6-2	M	6.6	1.2	1.4	10.0
SCH6-3	L	6.9	1.5	1.25	10.8
CH6S	Chin Style 6 Sizer Set				



Designed by Robert Gutstein M.D. Produces a cleft chin with proper dissection and suturing.

to order please call +91 91155 55555 or mail info@mitocon.in

HANSON
MEDICAL, INC.






CHIN STYLE 7

Cleft Chin

CAT #	Size	A	B	C	D
SCH7-1	S	6.0	1.2	0.7	7.2
SCH7-2	M	6.7	1.3	0.75	7.9
SCH7-3	L	7.0	1.5	0.95	8.3
CH7S	Chin Style 7 Sizer Set				

Designed by Stephen Giunta MD.

Ideal for creation of a feminine cleft chin. Requires proper dissection and suturing to produce a petite cleft.




CHIN STYLE 8

Cavity Back Chin

CAT #	Size	A	B	C	D
SCH8-1	S	6.0	1.6	0.6	9.7
SCH8-2	M	5.9	1.6	0.8	10.0
SCH8-3	L	5.9	1.7	0.9	9.7
SCH8-4	XL	5.9	1.8	1.0	10.0
CH8S	Chin Style 8 Sizer Set				

This chin provides a recessed area to accommodate the mons mentis prominence. Compare to *Implantech® FMG Series.




VOLOSHIN CHIN

Short Extended Anatomical

CAT #	Size	A	B	C	D
SVCH1	XS	5.4	1.4	0.4	6.9
SVCH2	S	5.4	1.4	0.6	7.1
SVCH3	M	5.4	1.4	0.8	7.1
SVCH4	L	5.9	1.55	1.0	7.5
SVCH5	XL	6.0	1.6	1.5	7.5
VCHS	Voloshin Chin Sizer Set				

Provides a recess for the mons mentis. It is designed by James Voloshin MD.





CHIN STYLE 10


Chin Jowl Style II

CAT #	Size	A	B	C	D
SCH10-1	S	5.9	1.4	0.6	10.3
SCH10-2	M	6.0	1.5	0.8	10.8
SCH10-3	L	6.0	1.5	0.9	10.9
SCH10-4	XL	6.0	1.5	1.1	11.0

CH10S Chin Style 10 Sizer Set

May be augmented with CH16

Provides lateral and forward projection. Compare to *Implantech® MPJC Series.



CHIN STYLE 11


Shelf Chin™

CAT #	Size	A	B	C	D
SCH11-1	S	7.4	1.30	0.55	10.2
SCH11-2	M	7.6	1.35	0.45	10.5
SCH11-3	L	7.8	1.40	0.55	10.8

CH11S Chin Style 11 Sizer Set

Increases height of chin

The only chin to offer a shelf to increase the vertical height of the chin as well as forward projection.




CHIN STYLE 12

Minimal Projection

CAT #	Size	A	B	C	D	E
SCH12-1	S	6.0	1.0	0.13	0.3	10.8
SCH12-2	M	6.0	1.0	0.13	0.4	10.9
SCH12-3	L	6.0	1.0	0.13	0.5	11.0
SCH12-4	XL	6.0	1.0	0.13	0.6	12.2

CH12S Chin Style 12 Sizer Set

Offers a scant 1mm forward projection, but augments the jowl groove. Compare to *Implantech® MPJ Series.






CHIN STYLE 13

Extended Anatomical

CAT #	Size	A	B	C	D
SCH13-1	S	6.35	1.05	0.30	11.8
SCH13-2	M	6.2	1.15	0.35	11.8
SCH13-3	L	6.25	1.25	0.50	12.5
SCH13-4	XL	6.1	1.3	0.55	12.5
CH13S	Chin Style 13 Size Set				




Provides anterior projection and extends tapering augmentation along the mandible toward the jaw angle. Compare to *Implantech® EAC Series.



CHIN STYLE 14

Magnum Chin Style I

CAT #	Size	A	B	C	D
SCH14-1	S	6.10	1.35	0.60	10.30
SCH14-2	M	6.40	1.50	0.66	11.00
SCH14-3	L	7.00	1.65	0.69	12.20
SCH14-4	XL	7.10	1.70	0.80	12.35
CH14S	Chin Style 14 Size Set				




Provides a square strong chin contour desired by male patients. Compare to *Implantech® TSCI Series.



CHIN STYLE 15

Magnum Chin Style II

CAT #	Size	A	B	C	D
SCH15-1	S	6.7	1.4	0.35	10.7
SCH15-2	M	8.2	1.6	0.40	12.7
SCH15-3	L	8.0	1.6	0.45	12.6
SCH15-4	XL	8.6	1.7	0.50	14.1
CH15S	Chin Style 15 Size Set				



Provides a greater square and masculine contour for the male patient. Compare to *Implantech® TSCII Series.



HANSON
MEDICAL, INC.

to order please call +91 91155 55555 or mail info@mitocon.in




CHIN STYLE 16


Add On Chin I

CAT #	A	B	C
SCH16-1	3.5	1.05	0.26
CH16S	Chin Style 16 Sizer Set		

For use with CH10



Choose this chin for increasing the projection of the SCH 10 Chin Jowl Style II.




CHIN STYLE 17

Add On Chin II

CAT #	A	B	C
SCH17-1	3.6	1.05	0.26
CH17S	Chin Style 17 Sizer Set		

For use with CH1 and CH13




Provides a 2mm increase in the projection of the SCH 1 and SCH 13



NASAL DORSAL

Lacks a Columella

CAT #	A	B	C
SND1-1F	4.25	1.0	0.15
SND1-2F	4.3	1.1	0.25
SND1-3F	4.4	1.1	0.35
SND1-4F	4.4	1.1	0.45




Eliminates the columella

to order please call +91 91155 55555 or mail info@mitocon.in

HANSON
MEDICAL, INC.






NASAL STYLES 1 THRU 16

CAT #	A	B	C	D	E
SN-1	5.6	2.3	0.85	0.65	0.6
SN-2	5.2	2.7	1.0	0.7	0.55
SN-3	4.5	2.7	0.9	0.6	0.45
SN-4	4.9	1.7	0.93	0.6	0.35
SN-5	4.9	1.9	1.05	0.75	0.45
SN-6	5.2	1.9	0.85	0.65	0.45
SN-7-1	5.5	2.5	0.95	0.5	0.35
SN-7-2	5.5	2.5	0.95	0.5	0.25
SN-8	4.7	2.5	1.05	0.65	0.35
SN-9	4.8	2.9	1.05	0.5	0.55


Hanson Medical offers the largest selection of nasal implant styles and sizes in the market. They are available in flesh tone or white. Most implants are easily modified with a scalpel. The columellar strut is almost always trimmed down, if not eliminated completely. Patient selection is critical. The patient should have adequate dermal and subcutaneous skin layers to cover the implant. Free sizers are available for each style of implant ordered.

CAT #	A	B	C	D	E
SN10-1	4.9	3.0	0.95	0.6	0.45
SN10-2	5.1	3.1	1.2	0.8	0.5
SN11	4.6	2.8	1.1	0.6	0.5
SN12-1	4.6	2.6	1.0	0.5	0.3
SN12-2	5.0	3.1	1.1	0.7	0.3
SN14	4.5	2.5	0.8	0.6	0.35
SN16	5.0	2.6	0.9	0.7	0.3




VOLOSHIN NASAL Dorsal Columella Design

CAT #	A	B	C	D	E
SV1	4.9	2.8	1.05	0.75	0.25
SV2	5.4	2.8	1.05	0.75	0.25
SV3	4.8	3.1	1.05	0.70	0.45
SV4	5.5	3.2	1.1	0.69	0.50
SV5	5.1	3.2	0.87	0.70	0.40




PECTORALIS 1**Medial Projection**


CAT #	Size	A	B	C
ACP1-1	S	14.0	9.0	1.6
ACP1-2	M	16.0	10.1	2.1
ACP1-3	L	16.8	11.4	2.3
ACP1-4	XL	18.0	14.4	2.9

**PECTORALIS 2****Anatomical Design**


CAT #	Size	A	B	C
ACP2-1	S	14.4	11.5	1.5
ACP2-2	M	15.2	12.5	2.0
ACP2-3	L	16.3	13.5	2.5
ACP2-4	XL	16.5	14.0	2.8

**PECTORALIS 3****Block Design**


CAT #	Size	A	B	C
ACP3-1	L	15.5	12.5	2.5

**GLUTEAL 4****Central Projection**


CAT #	Size	A	B	C
ACP4-1	M	14.8	10.3	2.7
ACP4-2	L	13.0	10.0	3.0

**GLUTEAL 5****Round**

CAT #	Size	A	B
ACP5-1	S	11.0	2.5
ACP5-2	L	13.0	3.0

**GLUTEAL 6 & 9****High Profile Projection**


CAT #	Size	A	B	C
ACP6-1	M-L	14.8	10.8	3.0
ACP9-1	M-L	14.0	11.0	3.5




We offer the Bircoll and Aiache designs in calf Implants. It is recommended that surgeons be conservative with the degree of augmentation and exercise caution concerning patient selection to avoid complications such as compartmental syndrome. Our Pectoral implants vary in projection widely and again we recommend a conservative approach in matching the size of the implant to the patient. The Style 2 pectoralis may be placed vertically or reversed and placed horizontally on the contralateral side. Therefore a right side can be used on the left side. One free size is available for each implant purchased.

CALF STYLES 1 & 2**Lateral & Medial Projections**


CAT #	Size	A	B	C
CCB1-1	XS	11.5	5.4	1.1
CCB1-2	S	12.7	6.1	1.3
CCB1-3	M	13.8	6.7	1.3
CCB1-4	L	14.4	7.5	1.3
CCB1-5	XL	16.2	8.5	1.5



CAT #	Size	A	B	C
CCB2-1	XS	11.5	4.6	2.0
CCB2-2	S	14.0	5.5	2.0
CCB2-3	M	18.0	6.0	2.0
CCB2-4	L	20.0	6.5	2.2
CCB2-5	XL	24.0	7.0	2.2



Hanson Medical is dedicated to providing you with High Quality Implants and outstanding customer service. Our customers are our first priority and we strive to accommodate your busy schedules. Most implants are in stock and ready to ship the same day.



All facial implants are provided STERILE.



All implant materials are FDA master filed solid silicone elastomer approved for long term implantation.

Facial implant dimensions are in centimeters and accurate to 1mm and approximately 30 durometer. Whereas, Pec, Glut and Calf implants have an approximate 7 durometer.

Size sets are available.

Custom Implants can be manufactured from a "Moulage Model" or by detailed drawings. Turn-around time is usually 7-10 days.

Moulage kits are available for Custom Implant sculpting.



INDIAN PARTNER :



MITOCCON BIOMED
INNOVATIONS GALORE

C/312, Lancelot Bldg, S. V. Road,
Borivali (W), Mumbai 400 092

M: +91 9892 444 444

E: info@mitocon.in W: www.mitocon.in

[Instagram](https://www.instagram.com/mitoconbiomed/) /info/inang [Facebook](https://www.facebook.com/mitoconbiomed/) /mitoconbiomed/ [YouTube](https://www.youtube.com/channel/UC...) /mitoconbiomed

HANSON
MEDICAL, INC.

PO Box 1160

Kingston, WA 98346

800-771-2215

360-297-1998 fax

www.hansonmedical.com