

CORPORATE OFFICE

A/204 - 211, THE ATLANTIC,
Opp. SRP Group IX, Makarpura Road,
Vadodara - 390 010, Guj. India


✉ info@orthomaxindia.net
admin@orthomaxindia.net

MFG. UNIT & REGD. OFFICE

C-1-B/886/4, G.I.D.C. Estate,
Makarpura, Vadodara - 390 010,
Guj. India

☎ +91 89800 15555
+91 89800 25555

WITH
30+
YEARS' EXPERIENCE



Elevate Your
Procedures to
A Whole New
Level



**CMF - PLASTIC
SURGERY
INSTRUMENTS**

www.orthomaxindia.net

Product Catalogue

www.orthomaxindia.net

PREFACE

It's a moment of great pleasure for us to present the all new edition of our products catalogue for Cranio-Maxillofacial (CMF) & Plastic surgery instruments.

This edition brings and demonstrates most comprehensive range of surgical instruments required in CMF & Plastic surgery.

We have been striving hard since 1993 to achieve excellence in manufacturing quality products at affordable prices for Orthopaedics, Cranio-Maxillofacial, Neuro & Plastic Surgery.

All our Instruments are approved by Indian FDA (CDSCO) under MD-5 License and manufactured under strict compliance of ISO 13485 standards which reflects our commitment and dedication towards the word "QUALITY".

The internationally acclaimed use of raw materials, standardized forging – machining- special polish treatment- cleaning & packaging in most hygienic condition enrich our products to "PERFECTION". We keep proper check and measure every little details count in making every product fit for critical surgical procedure. We aim at making our products most vital part of your surgery because we believe that the quality of the products elevates your procedures to a whole new level. We thank all our users and channel partners for their valuable feedback and suggestions for the continual improvement to achieve the highest parameters of product quality & safety.

All our products can be easily bought by contacting our nearest Dealers or by simply sending your purchase enquiry to customercare@orthomaxindia.net or Whatsapp to: +91 89800 25555

For further details, please spare time to visit www.orthomaxindia.net



2024 EDITION


CONTENTS

BASIC BONE PLATING INSTRUMENTS	1-5
BONE HOLDING FORCEPS	6-7
RETRACTORS	7-12
ELEVATORS	12-14
BONE AWLS & HOOKS	15
BONE CUTTERS & NIBBLERS	16
OSTEOTOMES-CHISELS-GOUGES-CURRETTES	17-18
MOUTH GAGS	19
CLEFT LIP & PALATE INSTRUMENTS	19-20
EXTERNAL DISTRACTORS & FIXATORS	21
BOXES & CONTAINERS	22
OTHER INSTRUMENTS	23-25
GENERAL INSTRUMENTS - PREMIUM	26-28
BASIC TRAUMA KIT	29
IMPACTION KIT	30
INSTRUMENTS SET FOR ORTHOGNATHIC / MAXILLOFACIAL SURGERY	31
INSTRUMENTS SET FOR CLEFT LIP & PALATE	32
INSTRUMENTS SET FOR RHINOPLASTY	33
INSTRUMENTS SET FOR OSTEOTOMY & BONE GRAFTING	34
INSTRUMENTS SET FOR WIRING	35
INSTRUMENTS SET FOR MICRO SURGERY	36
BASIC SET FOR PLASTIC SURGERY	37
MINOR ORAL SURGERY KIT	38
INSTRUCTIONS FOR USE & MAINTENANCE OF INSTRUMENTS	39

GUARANTEES

Ortho Max Guarantees the quality for all its instruments. Any item which is found to be defective because of material or poor craftsmanship will be repaired or replaced free of cost within six months from date of supply.

CDSO LICENCE TO MANUFACTURE CLASS A & B DEVICES



www.orthomaxindia.net

BASIC BONE PLATING INSTRUMENTS

Bone Elevator - Hohman Type



Ref.	Size
252.006	6mm
252.008	8mm

Plier



Ref.
253.001

Depth Gauge



Ref.	Size
254.150	1.5mm
254.200	2.0mm
254.250	2.5mm

Trocar Sleeve / Drill Guide



Ref.	Size
255.001	2.35mm

Self Holding Screw Driver

For Single Slotted Screws



Ref.	Size
257.150	1.5mm
257.200	2.0mm
257.250	2.5mm

* Available in Aluminum handle also

Ordinary Screw Driver

For Single Slotted Screws



Ref.	Size
258.150	1.5mm
258.200	2.0mm
258.250	2.5mm

Cross Drive® Screw Driver

Fixed Handle



Ref.	Size
493.120	1.2 mm (Non Locking)
493.150	1.5 mm (Non Locking)
493.180	1.8 / 2 mm (Locking)
493.200	2.0 mm (Non Locking)
493.230	2.3 / 2.5 mm (Locking)
493.250	2.5 mm (Non Locking)

Universal Cross Drive® Screw Driver



Ref.
518.001

Cross Drive® Screw Driver Handle
Quick Coupling



Ref.
493.001

Ultralock® Screw Driver Shaft



Ref.	Size
493.118	1.8 mm / 2mm
493.123	2.3 mm / 2.5 mm

Cross Drive® Screw Driver Shaft



Ref.	Size
493.012	1.2mm



Ref.	Size
493.015	1.5mm



Ref.	Size
493.020	2.0mm



Ref.	Size
493.025	2.5mm

Bone Tap - S.S.



Ref.	Size
259.150	1.5 mm
259.200	2.0 mm
259.250	2.5 mm

iDrill® Screw Driver Handle



Ref.
1070.001

iDrill® Screw Driver Shaft




Ref.	Size
1071.001	1.5 /2.0mm

iDrill® Screw Driver Shaft with Sleeve



Ref.	Size
1071.002	1.5 /2.0mm

Mini Bone Drill



Ref.
260.002

Plate Bending Forceps
Modified



Ref.
261.001

Plate Bending Forceps with Pins
For without bar plates



Ref.	Size
261.150	1.5 mm
261.200	2.0 mm
261.250	2.5 mm

Modelling Lever



Ref.
262.001

Screw / Plate Holding Forceps



Ref.
521.001

Plate Holding Fork



Ref.
263.002

Counter Sink for Lag Screw



Ref.	Size
564.200	2.0 mm
564.250	2.5 mm

Template for Ultra Lock® Reconstruction Plate
Primary



Straight	Ref.	Size
	463.116	16 Holes
	463.120	20 Holes
	463.124	24 Holes

Single Angle	Ref.	Size
	463.220	20 Holes
	463.224	24 Holes

Double Angle	Ref.	Size
	463.328	4+20+4 Holes
	463.334	6+22+6 Holes

Template for *Ultra Lock*® Reconstruction Plate
Regular & Secondary



Straight	Ref.	Size
	463.416	16 Holes
	463.420	20 Holes
	463.424	24 Holes

Single Angle	Ref.	Size
	463.520	20 Holes
	463.524	24 Holes

Double Angle	Ref.	Size
	463.628	4+20+4 Holes
	463.634	6+22+6 Holes

Drill Sleeve for *Ultra Lock*® Reconstruction Plate



Ref.	Size
490.120	1.2mm
490.150	1.5mm
490.170	1.7mm
490.200	2mm

Mandibular Reconstruction Plate Bender
Roller Type



Ref.
265.001

Plate Bender in Pair -Small
for *Ultra Lock*® reconstruction plate



Ref.
571.003

Multi Plate Cutter
Circular Type



Ref.
506.001

Reconstruction Multi Angle Plate Bender



Ref.
526.001

Arch Bar Cutter




Ref.
451.001

Wire Cutter
Ordinary



Ref.
450.004

Wire / Plate Cutter With TC Tip



Ref.	Size
549.002	Small
549.001	Big

Wire Twisting Forceps



Ref.
446.001



Ref.
446.002

Mini Plate Cutter
Vertical



Ref.
452.002

Wire Twister
Pencil Type



Ref.
447.001

Mesh Cutter



Ref.
520.001

Mesh Bender



Ref.
454.001

Plate Cutter
Mini



Ref.
452.001

Transbuccal Device with Trocar & Cheek Retractor



Ref.
438.001

BONE HOLDING FORCEPS


www.orthomaxindia.net

Rowe Maxillary Disimpaction Forceps



Ref.	Size
266.001	Right & Left in Pair

Hayton William Forceps



Ref.	Size
268.01F	Forward Traction
268.01D	Downward Traction

Asch Nasal Septum Forceps



Ref.	Size
269.01S	Small



Ref.	Size
269.01B	Big

Walsham Nasal Septum Forceps

Right & Left in Pair



Ref.
270.001

Reduction Bone Holding Forceps



Pointed

Ref.
273.006



Slightly Curved

Ref.
273.003

www.orthomaxindia.net

Reduction Bone Holding Forceps

Serrated with Ratchet



Ref.
275.001

TMJ Spreading Forceps



Ref.
277.001

Bone & Plate Holding Forceps - Guarded Point

Right & Left in Pair



Ref.
274.001

Fibula Holding Forceps



Ref.
916.001

RETRACTORS

Czerny Retractor



Ref.
281.001

Double Hook Skin Retractor



Ref.
282.001

McIndoe Retractor



Ref.
283.001

Rayne Mandibular Retractor



Ref.
285.001

Davis Double Ended Retractor



Ref.
286.001

Forked Ramus Retractor



Ref.
288.001

Kilner Retractor Double Ended



Ref.	Size
290.001	Small
290.002	Medium
290.003	Big

Kilner Nasal Retractor



Ref.
292.001

Aufricht Nasal Retractor



Solid

Ref.
293.001



Fenestrated

Ref.
293.002

Chin Retractor



Ref.
287.001

Channel Retractor



Ref.	Size
289.008	8mm
289.012	12mm

Channel Retractor

Big Handle



Ref.
289.0XB

Kilner Skin Retractor

Double Ended



Ref.
294.001

Aufricht Nasal & Double Hook Retractor



Ref.
498.001

Malleable Retractor - Copper




Ref.	Size
296.008	8" X 10mm
296.010	10" X 15mm
296.012	12" X 20mm

Desmarres Retractor



Ref.	Size
301.001	Small
301.002	Medium
301.003	Big

Austin Tissue Retractor



Regular

Ref.
298.001



Half Moon

Ref.
298.002



Forked

Ref.
298.003

Self Retaining Skin Retractor - "Alm's"



Small

Ref.
299.001



Big

Ref.
299.002

Right Angle Retractor

Lengenback



Ref.	Size
300.620	Small 6 X 20 mm



Ref.	Size
300.1030	Med. 10 X 30 mm



Ref.	Size
300.1545	Big 15 X 45 mm

Condyle Retractor Double Ended



Ref.
303.002

Kilner Alae Retractor




Ref.
307.001

Lip & Cheek Retractor Bishop



Ref.
305.001

Self Retaining Mastoid Retractor



Ref.	Size
309.004	4"
309.005	5"
309.006	6"

Volkman Retractor



Ref.
308.001

Browne Lingual Flap Retractor



Ref.
310.001

TMJ Retractor



Ref.	Size
311.610	6mm X 10mm

Condylar Neck Retractor Right & Left in Pair



Ref.
304.001

Condylar Alveolar Retractor



Ref.
312.001

Rowe Orbital Floor Retractor Right & Left in Pair



Ref.
313.001

Soft Tissue Retractor Reverse Lengenback



Ref.
315.001

Soft Tissue Retractor Concave Blade



Ref.
316.002

SENN Retractor



Ref.
318.001

Minnesota Retractor



Ref.
317.001

Superior Ramus Separator Smith



Ref.
314.001

Universal Maxillary Spreader



Ref.
314.002

Sigmoid Notch Retractor



Ref.
321.001

Fibula Retractor



Ref.
524.001

RETRACTORS

www.orthomaxindia.net

Intraoral Retractor

3rd Molar



Ref.
323.001

Plastic Cheek Retractor



Ref.	Size
324.001	Small
324.002	Medium
324.003	Big

Ramus Retractor

Big Handle



Ref.	Size
328.020	20mm
328.025	25mm
328.030	30mm

Deaver Retractor



Ref.	Size
494.001	1/2"
494.002	1"

ELEVATORS

Mcdonald Dissector



Ref.
330.001

Apexo Elevator

Right & Left in Pair



Ref.
331.101

Ward Periosteal Elevator

Fan Shaped



Ref.
332.001

Dental Elevator

Straight



Ref.
331.001

Right & Left in Pair



Ref.
331.002

Bristow Elevator



Ref.
334.001

www.orthomaxindia.net

Periosteal Elevator

Fibre Handle



Ref.
333.001



Ref.
333.003



Ref.
333.005



Ref.
333.007



Ref.
333.002



Ref.
333.004



Ref.
333.006



Ref.
333.008

Cryer Elevator

Right & Left in Pair



Ref.
335.001

Cryer X Bar Elevator

Right & Left in Pair



Ref.
336.001

Warwick James (Ward) Elevator

Straight



Ref.
337.001

Right & Left in Pair



Ref.
337.002

Rowe Zygomatic Elevator



Ref.	Size
338.008	8mm
338.012	12mm

Kilner Zygomatic Elevator



Ref.
339.001

Farabeuf Rugine



Ref.	Size
340.001	Straight
340.002	Curved

Periosteal Elevator

Double Ended



Ref.
342.001

Septum Elevator



Ref.
345.001

Hajeck Septum Elevator



Ref.
349.001

Molt Elevator



Ref.	Size
351.001	Big
351.002	Small

Cottle Elevator

Graduated



Ref.
352.001

Howarth Elevator



Ref.
341.001

Dissector Probe

Double Ended



Ref.
344.001

Alveolar Dissector

Double Ended



Ref.
348.001

Dingman Elevator



Ref.
1067.001

Coupland Elevator



Ref.	Size
350.003	3 mm
350.004	4 mm

Root Tip Elevator



Ref.
353.001

Kelsey Fry Bone Awl



Ref.
356.001



Ref.
356.002

Zygomatic Bone Awl



Ref.	Size
354.006	6"
354.008	8"

Gillies Skin Hook



Ref.
358.001

Kilner Skin Hook




Ref.
360.001

Poswillo Malar Hook



Ref.
362.001

Mandibular Awl Set of 5



Ref.
355.001

Flat Handle



Ref.
358.001

Round Handle



Ref.
358.002

Double (McIndoe) Skin Hook



Ref.
361.001

Stacey Zygomatic Hook



Ref.
364.001

BONE CUTTERS & NIBBLERS

www.orthomaxindia.net

Bone Cutting Forceps

Single Action



Ref.	Size
365.S7S	7.5" Straight
365.S7C	7.5" Curved

Bone Cutting Forceps

Double Action



Ref.	Size
365.D5S	5.5" Straight
365.D5C	5.5" Curved
365.D7S	7.5" Straight
365.D7C	7.5" Curved
365.D9S	9" Straight
365.D9C	9" Curved

Strong Bone Cutter



Ref.
367.001

Rib Cutter

Big



Ref.
369.001

Bone Nibbler

Single Action



Ref.	Size
366.S5S	5.5" Straight
366.S5C	5.5" Curved
366.S7S	7.5" Straight
366.S7C	7.5" Curved

Bone Nibbler

Double Action



Ref.	Size
366.D5S	366.D5S
366.D5C	366.D5C
366.D7S	366.D7S
366.D7C	366.D7C
366.D9S	366.D9S
366.D9C	366.D9C


Doyen's Rib Shear



Ref.
368.001

Rib Raspatory

Right & Left in Pair



Ref.
370.001

www.orthomaxindia.net

OSTEOTOMES-CHISELS-GOUGES-CURRETTES

Fine Gouge



Ref.	Size
371.002	2mm
371.003	3mm
371.004	4mm
371.006	6mm
371.008	8mm
371.010	10mm

Fine Osteotome



Straight

Ref.	Size
373.S02	2mm
373.S03	3mm
373.S04	4mm
373.S05	5mm
373.S06	6mm
373.S08	8mm
373.S10	10mm
373.S12	12mm



Curved

Ref.	Size
373.C02	2mm
373.C03	3mm
373.C04	4mm
373.C05	5mm
373.C06	6mm
373.C08	8mm
373.C10	10mm
373.C12	12mm

Cottle Osteotome

Half Moon



Ref.	Size
375.001	8mm
375.002	10mm

Tessier Osteotome



Straight

Ref.	Size
372.S05	5mm
372.S10	10mm
372.S15	15mm
372.S20	20mm



Curved

Ref.	Size
372.C05	5mm
372.C10	10mm
372.C15	15mm
372.C20	20mm

Fine Chisels



Ref.	Size
374.002	2mm
374.003	3mm
374.004	4mm
374.005	5mm
374.006	6mm
374.008	8mm
374.010	10mm
374.012	12mm

Nasal Septum Osteotome with Guard

"U" Shape



Ref.
377.001

Nasal Chisel with Guard

Right & Left in Pair



Ref.	Size
378.004	4mm
378.005	5mm
378.006	6mm
378.008	8mm

Pterigoidous Chisel



Ref.	Size
380.006	6mm
380.008	8mm
380.010	10mm
380.012	12mm

Volkman Bone Scoop

Double Ended



Ref.
384.002

Bone Scoop

Angled



Ref.
387.001

Ramus Pusher



Ref.
379.002

Ramus Stripper



Ref.
379.008

Bone Grafting Gouge



Ref.	Size
383.002	2mm
383.004	4mm
383.006	6mm
383.008	8mm
383.010	10mm

Lucas Currettes

Double Ended



Ref.
386.001

Heister Jaw Opener



Ref.	Size
393.001	Adult
393.003	Paediatric

Doyen Mouth Gag



Ref.	Size
390.001	Adult
390.002	Paediatric

Cleft Palate Raspatory - French Pattern

Right & Left in Pair



Ref.
395.001

McIndoe Cleft Palate Raspatory



Ref.	Size
398.001	Small
398.002	Medium
398.003	Big

Ferguson Mouth Gag



Ref.
394.001

Dingman Mouth Gag with 3 Blades



Ref.	Size
391.001	Adult
391.002	Paediatric

CLEFT LIP & PALATE INSTRUMENTS

Cleft Palate Raspatory - Barsky Type

Right & Left in Pair



Ref.
396.001



Ref.
396.002

Cleft Palate Raspatory

Curved



Ref.	Size
397.001	Up
397.002	Down

Barsky Pharyngeal Flap Elevator



Ref.
399.001


Cleft Palate Knife

Triangular



Ref.
401.001

Cronin Cleft Palate Elevator



Ref.	Size
403.01R	Small - Right
403.01L	Small - Left
403.02R	Med - Right
403.02L	Med - Left
403.03R	Big - Right
403.03L	Big - Left

Spoon Shape Elevator

Double Ended



Ref.
405.001

Elevator Double Ended



Ref.	Size
407.001	Right & Left
407.002	Side & Flat

Mitchell Trimmer



Ref.
400.001

Cleft Palate Hook



Right Angle

Ref.
402.001


Down Bent



Ref.
402.002

Palate Dissector

Right & Left in Pair



Ref.
404.001

Miller Bone File



Ref.	Size
406.001	Straight
406.002	Curved

External Mandible Distractor

Uni-Directional



Ref.
415.001

Intended Use	For Distraction Osteogenesis of Mandible Body
Rotation	Anti-Clock wise two full turns for 1mm distraction
Pin Size	2.35mm (St. Hand Piece Fit)

External Fixator For Finger



Ref.	Size
416.001	1.5mm X 3"
416.002	2mm X 4"

Activator for External Mandible Distractor

Box Type



Ref.
417.004

Intended Use	For Rotation of External Mandible Distractor
--------------	--

Activator for Intraoral Distractor

Rigid



Ref.
516.001

Intended Use	For Rotation of Intraoral Distractor
--------------	--------------------------------------

External Mandible Distractor

Bi-Directional



Ref.
415.002

Intended Use	For Distraction Osteogenesis of Mandible Body
Rotation	Anti-Clock wise two full turns for 1mm distraction
Pin Size	2.35mm (St. Hand Piece Fit)

External Mandible Distractor Guide Sleeve

2.35mm



Ref.
542.001

Intended Use	For Parallel Pin Insertion
--------------	----------------------------

Wrench



Ref.
791.001

Intended Use	For Stopper nut fitting of External Mandible Distractor
--------------	---

Screw Driver for External Mandible Distractor

Hexagonal



Ref.
417.002

Intended Use	For Schanz Screw Locking
--------------	--------------------------

Activator for Intraoral Distractor

Hinged



Ref.
517.001

Intended Use	For Rotation of Intraoral Distractor
--------------	--------------------------------------

BOXES & CONTAINERS

www.orthomaxindia.net

Skin Grafting Handle Container



Ref.
418.002

Drill Bits Container



Ref.	Size
424.001	4"

Mini Containers for Implants



Ref.
421.150



Ref.
421.200



Ref.
421.250

General Instruments Box with Silicon Mat



Ref.	Size
428.002	L 24cm x W 20cm x H 2.5cm

K Wire Container

Aluminum



Ref.	Size
425.206	6"

K Wire Container Aluminum with Partition



Ref.	Size
425.306	6"

Aluminum General Box

With Tray



Small

Ref.	Size
492.001	L 24cm x W 16.5cm x H 7cm



Medium

Ref.	Size
492.002	L 13cm x W 13cm x H 10cm



Big

Ref.	Size
492.003	L 37cm x W 16cm x H 13cm

www.orthomaxindia.net

Bone File/Rasp



Ref.
430.001

Tendon Tunneling Forceps



Ref.	Size
432.08S	8" Straight
432.08C	8" Curved
432.10S	10" Straight
432.10C	10" Curved

Cartilage Crusher



Rectangular

Ref.
434.001

Round



Ref.
434.002

Gigli Saw Handle

in Pair



Ref.
436.001

Tessier Maxillary Mobiliser

Right & Left in Pair



Ref.
439.001

OTHER INSTRUMENTS

Killien Nasal Speculum



Ref.
431.001

Tendon Stripper



Ref.	Size
433.416	4mm x 16"
433.516	5mm x 16"
433.616	6mm x 16"

Gigli Saw Wire



Ref.
435.001

Gigli Saw Guide



Ref.
437.001

Nasal Saw



Straight

Ref.
440.001



Left

Ref.
440.002



Right

Ref.
440.003

Finger Saw



Ref.
441.001

OTHER INSTRUMENTS

www.orthomaxindia.net

Fickling Dental Mallet



Ref.
443.001

Nylon Faced Hammer



Ref.	Size
444.001	Small
444.002	Medium
444.003	Big

Caliper Castroviejo



20mm



40mm

Ref.	Size
448.001	Short-Straight (20mm)
448.002	Short-Curved (20mm)
448.003	Long-Straight (20mm)
448.004	Long-Curved (20mm)
448.005	Long-Straight (40mm)
448.006	Long-Curved (40mm)

K Wire Bender/Cutter/Plier - 3 in 1



Ref.	Size	Cutting Capacity
453.001	Small	Up to 1.5mm Dia
453.002	Big	Up to 2mm Dia

Plier (for EzyBar)




Ref.
253.002

Plate Spacer (for EzyBar)



Ref.
527.001

K Wire Introducer



Ref.
445.001

Suction Cannula

Mastoid



Ref.	Size	Ref.	Size
455.100	No. 1	455.400	No. 4
455.200	No. 2	455.500	No. 5
455.300	No. 3	455.600	No. 6

Suction Cannula

Frazier



Ref.	Size	Ref.	Size
489.100	No. 1	489.400	No. 4
489.200	No. 2	489.500	No. 5
489.300	No. 3	489.600	No. 6

www.orthomaxindia.net

Suction Cannula

Frazier with Ball Tip



Ref.	Size
508.100	No. 1
508.200	No. 2
508.300	No. 3
508.400	No. 4
508.500	No. 5
508.600	No. 6

Bone Trephine




Ref.	Size
458.002	4.5 mm Graft

Skin Grafting Handle



Ref.	Size
487.0BR	Big - Right
487.0BL	Big - Left

Tongue Depressor



Fenestrated

Serrated

Ref.
457.001

Ref.
457.002

Mouth Props



Ref.	Size
460.001	Set of 3
460.002	Set of 4

Carroll Girard Screw



Ref.
1063.001

Scale Simple

S.S.



Ref.	Size
563.006	6"

Drill/Screw/Wire Gauge Scale



Ref.	Size
459.001	6"

Scissors - Fine Irish Type



Ref.	Size
477.I3S	3.5" - Straight
477.I3C	3.5" - Curved
477.I4S	4.5" - Straight
477.I4C	4.5" - Curved

Scissors - Metzenbaum



Ref.	Size
477.M4S	4" - Straight
477.M4C	4" - Curved
477.M5S	5" - Straight
477.M5C	5" - Curved
477.M6S	6" - Straight
477.M6C	6" - Curved
477.M7S	7" - Straight
477.M7C	7" - Curved

Scissors - "Heath" Stitch Cutting



Ref.	Size
477.H01	4"
477.H02	6"

Scissors - Mayo



Ref.	Size
477.Y6S	6.5" - Straight
477.Y6C	6.5" - Curved
477.Y7S	7.5" - Straight
477.Y7C	7.5" - Curved
477.Y8S	8.5" - Straight
477.Y8C	8.5" - Curved

Scissors - Stevens Tenotomy



Ref.	Size
477.T4S	4" - Straight
477.T4C	4" - Curved
477.T5S	5" - Straight
477.T5C	5" - Curved
477.T6S	6" - Straight
477.T6C	6" - Curved
477.T7S	7" - Straight
477.T7C	7" - Curved

Scissors - "Northbent" Suture Cutting



Ref.
477.S01

Scissors - Kilner



Ref.
477.K05

Scissors - Angular for Cleft Palate



Ref.
477.A01

Scissors - Micro / Castroviejo



Ref.	Size
477.C6S	6" - Straight
477.C6C	6" - Curved

Needle Holder - Ryder

For 5.0/6.0



Ref.	Size
478.305	5"
478.306	6"
478.308	8"

Needle Holder - Deep Serration

For 3.0/4.0



Ref.	Size
478.005	5"
478.006	6"
478.008	8"

Needle Holder

Bozeman Angular



Ref.	Size
478.106	6"
478.108	8"

Micro Needle Holder

Castroviejo



Ref.	Size
478.206S	6" - Straight
478.206C	6" - Curved

Crile Artery Forceps



Ref.	Size	Ref.	Size
479.06S	6" - Straight	479.08S	8" - Straight
479.06C	6" - Curved	479.08C	8" - Curved
479.07S	7" - Straight	479.10S	10" - Straight
479.07C	7" - Curved	479.10C	10" - Curved

Mixer Artery Forceps




Ref.	Size
479.16C	6"
479.18C	8"

Mosquito Artery Forceps



Ref.	Size
480.04S	4" - Straight
480.04C	4" - Curved
480.05S	5" - Straight
480.05C	5" - Curved

Dandy Artery Forceps



Ref.	Size
522.005	5.5"

Dissecting Forceps

Adson



Ref.	Size
481.05P	5" - Plain
481.05T	5" - Tooth
481.07P	7" - Plain
481.07T	7" - Tooth

Dissecting - Micro Forceps



Ref.	Size
481.16P	6" - Plain
481.26T	6" - Tooth

BabCock Forceps



Ref.	Size
528.006	6"
528.008	8"

Vessel Dilators



Ref.	Size
525.001	5"

Allis Tissue Holding Forceps



Ref.	Size
482.004	4"
482.006	6"
482.008	8"
482.010	10"

Kocher Clamps



Ref.	Size	Ref.	Size
484.06S	6" - Straight	484.08S	8" - Straight
484.06C	6" - Curved	484.08C	8" - Curved

Dissecting Forceps

Adson Tooth - Angled



Ref.
1068.06T

Jeweller's Forceps



Ref.	Size
519.001	5" - Plain
519.002	5" - Tooth

Sponge Holding Forceps



Ref.	Size
483.006	6"
483.008	8"
483.010	10"

Towel Forceps

Guarded



Ref.	Size
485.005	5.5"

B. P. Handle

Single Ended



Flat Handle

Ref.	Size
486.003	No. 3
486.004	No. 4
486.007	No. 7



Round Handle









Ref.	Size
486.103	No. 3
486.104	No. 4

B. P. Handle

Double Ended



Ref.	Size
486.034	3 x 4

 Ref. 356.002	 Ref. 262.001	 Ref. 258.150	 Ref. 258.250
 Ref. 258.200	 Ref. 253.001	 Ref. 273.006	 Ref. 411.100
 Ref. 411.150	 Ref. 411.200	 Ref. 563.006	 Ref. 257.150
 Ref. 257.250	 Ref. 257.200	 Ref. 450.004	 Ref. 446.001




Ref. 492.001

Ref.	Instrument Contents	Qty
492.001	Container For Basic Trauma Kit	01 No
356.002	Kelsey Fry Bone Awl - Fibre Handle	01 No
262.001	Modelling Lever	01 No
258.150	Ordinary Screw Driver F/H 1.5mm	01 No
258.200	Ordinary Screw Driver F/H 2mm	01 No
258.250	Ordinary Screw Driver F/H 2.5mm	01 No
253.001	Plier	01 No
273.006	Reduction Bone Holding Forceps Pointed 6"	01 No
411.100	S.S. Drill Bits For M. M. Regular 1 x 12 x 65mm	01 No

Ref.	Instrument Contents	Qty
411.150	S.S. Drill Bits For M. M. Regular 1.5 x 16 x 65mm	01 No
411.200	S.S. Drill Bits For M. M. Regular 2 x 20 x 65mm	01 No
563.006	Scale Simple - S.S. 6"	01 No
257.150	Self Holding Screw Driver F/H 1.5mm	01 No
257.200	Self Holding Screw Driver F/H 2mm	01 No
257.250	Self Holding Screw Driver F/H 2.5mm	01 No
450.004	Wire Cutter - Ordinary	01 No
446.001	Wire Twisting Forceps (Thick)	02 Nos


IMPACTION KIT

www.orthomaxindia.net



Ref.	Instrument Contents	Qty
486.003	B. P. Handle No. 3	01 No
341.001	Howarth Elevator	01 No
337.001	Warwick James Elevator Straight	01 No
337.002	Warwick James Elevator - In pair Left & Right	01 Pair
351.001	Molt Elevator	01 No
298.001	Austin Tissue Retractor - Regular	01 No
350.004	Coupland Elevator - 4mm	01 No
353.001	Root Tip Elevator	01 No
386.001	Lucas Currette D/E	01 No
460.001	Mouth Props Set of 3	01 Set
1052.E106	Eco Needle Holder Deep Serration 6"	01 No
477.S01	Scissors - "Northbent" Suture Cutting	01 No
406.001	Miller Bone File	01 No
489.400	Suction Cannula - Frazier No. 4	01 No
373.S04	Fine Osteotome Straight 4 mm	01 No
1055.E6T	Eco Dissecting Forceps - Adson-6" Tooth	01 No
1054.E05C	Eco Mosquito Artery Forceps - 5" Curved	01 No
492.001	Container for Impaction Kit	01 No

INSTRUMENTS SET FOR ORTHOGNATHIC / MAXILLOFACIAL SURGERY




Ref.	Instrument Contents	Qty
266.001	Rowe Maxillary Disimpaction Forceps - Right & Left	01 Pair
268.01F	Hayton William Forceps - Forward Traction	01 No
273.006	Reduction Bone Holding Forceps Pointed 6"	01 No
277.001	TMJ Spreading Forceps	01 No
282.001	Double Hook Skin Retractor	01 No
285.001	Rayne Mandibular Retractor	01 No
287.001	Chin Retractor	01 No
288.001	Forked Ramus Retractor	01 No
289.012	Channel Retractor - 12mm	01 No
290.003	Kilner Retractor D/E "C" Shape - Big	01 No
298.001	Austin Tissue Retractor - Regular	01 No
300.620	Right Angle Retractor Small 6 x 20 mm	01 No

Ref.	Instrument Contents	Qty
300.1030	Right Angle Retractor Medium 10 x 30 mm	01 No
300.1545	Right Angle Retractor Big 15 x 45 mm	01 No
303.002	Condyle Retractor Double Ended	01 No
309.006	Self Retaining Mastoid Retractor 6"	01 No
313.001	Rowe Orbital Floor Retractor - Right & Left	01 Pair
314.001	Superior Ramus Separater (Smith)	01 No
315.001	Soft Tissue Retractor - Reverse Lengenback	01 No
334.001	Bristow Elevator	01 No
333.003/5	Periosteal Elevator Fibre Handle 2 Shapes	02 Nos
338.012	Rowe Zygomatic Elevator 12mm	01 No
337.002	Warwick James Elevator - Right & Left	01 Pair
492.003	Aluminium General Box - Big	01 No

INSTRUMENTS SET FOR CLEFT LIP & PALATE

www.orthomaxindia.net




Ref.	Instrument Contents	Qty
391.001	Dingman Mouth Gag with 3 Blades Adults	01 Set
395.001	Cleft Palate Raspatory - French Pattern Right & Left	01 Pair
396.002	Cleft Palate Raspatory - Barsky Type Double Eneded	01 No
397.001	Cleft Palate Raspatory - Curved Up	01 No
397.002	Cleft Palate Raspatory - Curved Down	01 No
398.002	McIndoe Cleft Palate Raspatory - Medium	01 No
400.001	Mitchell Trimmer	01 No
401.001	Cleft Palate Knife Triangular	01 No
402.001-2	Cleft Palate Hook Right Angle & Down Bent	02 Nos
403.02R	Cronin Cleft Palate Elevator Medium Right	01 No
404.001	Palate Dissector - Right & Left	01 Pair
405.001	Spoon Shape Elevator - Double Ended	01 No
406.002	Miller Bone File - Curved	01 No
407.001	Elevator Double Ended - Angled Right & Left	01 No
407.002	Elevator Double Ended - Angled Side & Flat	01 No
341.001	Howarth Elevator	01 No
294.001	Kilner Skin Retractor - Double Ended	01 No
492.001	Container for Cleft Lip & Palate Set	01 No



INSTRUMENTS SET FOR RHINOPLASTY

www.orthomaxindia.net




Ref.	Instrument Contents	Qty
358.001	Gillies Skin Hook Flat Handle	01 No
358.002	Gillies Skin Hook Round Handle	01 No
360.001	Kilner Skin Hook	01 No
282.001	Double Hook Skin Retractor	01 No
352.001	Cottle Elevator - Graduated	01 No
340.002	Farabeuf Rugins Curved	01 No
345.001	Septum Elevator	01 No
293.001	Aufricht Nasal Retractor - Solid	01 No
307.001	Kilner Alae Retractor	01 No
318.001	SENN Retractor	01 No
286.001	Davis Double Ended Retractor	01 No
374.002-4	Fine Chisels 2, 4mm	02 Nos

Ref.	Instrument Contents	Qty
377.001	Nasal Septum Osteotomes With Guard 'U' Shape	01 No
375.001	Cottle Osteotome (Halfmoon) - 8mm	01 No
373.S04-06	Fine Osteotome 4, 6mm - St.	02 Nos
378.004	Nasal Chisel With Guard 4mm - Right & Left	01 Pair
371.003	Fine Gouge 3mm	01 No
440.001	Nasal Saw - St	01 No
349.001	Hajeck Septum Elevator	01 No
430.001	Bone File Double Ended	01 No
269.01B	Asch Nasal Septum Forceps Big	01 No
270.001	Walsham Nasal Septum Forceps - Right & Left	01 Pair
492.002	Aluminium General Box - Medium	01 No

INSTRUMENTS SET FOR OSTEOTOMY & BONE GRAFTING

www.orthomaxindia.net



Ref.	Instrument Contents	Qty
333.005	Periosteal Elevator Fibre Handle (No. 5)	01 No
365.D5C	Bone Cutting Forceps Double Action Curved 5.5"	01 No
366.D5C	Bone Nibbler Doble Action Curved 5.5"	01 No
372.S05-10-15	Tessier Osteotome St. 5, 10, 15 mm	03 Nos
372.C10	Tessier Osteotome Cd. 10 mm	01 No
373.S04-06-08	Fine Osteotomes St. 4, 6, 8 mm	03 Nos
373.C06	Fine Osteotomes Cd. 6 mm	01 No
374.006	Fine Chisels 6mm	01 No
444.002	Nylon Faced Hammer Medium	01 No
448.001	Calipers Castroviejo - Short 20mm St.	01 No
458.002	Bone Trephine 4.5 mm	01 No
434.001	Cartilage Crusher Rectangular Shape	01 No
380.010	Pterigoidous Chisel 10 mm	01 No
492.002	Aluminium General Box - Medium	01 No

www.orthomaxindia.net

INSTRUMENTS SET FOR WIRING



Ref.	Instrument Contents	Qty
354.006	Zygomatic Bone Awl 6"	01 No
356.002	Kelsey Fry Bone Awl Fibre Handle	01 No
355.001	Mandibular Awl Set of 5	01 Set
362.001	Poswillo Malar Hook	01 No
364.001	Stacey Zygomatic Bone Hook	01 No

Ref.	Instrument Contents	Qty
446.001	Wire Twisting Forceps - Thick	01 No
450.001	Wire Cutter Regular	01 No
451.001	Arch Bar Cutter	01 No
492.001	Container for Wiring Set	01 No

INSTRUMENTS SET FOR MICRO SURGERY

www.orthomaxindia.net




Ref.	Instrument Contents	Qty
478.206S	Needle Holder Castroviejo - 6" Straight	01 No
478.206C	Needle Holder Castroviejo - 6" Curved	01 No
481.16P	Micro Forceps 6" Plain	01 No
481.26T	Micro Forceps 6" Tooth	01 No
477.C6S	Scissors - Micro / Castroviejo Straight	01 No

Ref.	Instrument Contents	Qty
477.C6C	Scissors - Micro / Castroviejo Curved	01 No
519.001	Jeweller's Forceps	01 No
525.001	Vessel Dilator	01 No
428.002	General Instruments Box with Silicone Mat	01 No

BASIC SET FOR PLASTIC SURGERY

www.orthomaxindia.net



Ref.	Instrument Contents	Qty
358.002	Gillies Skin Hook - Round Handle	02 Nos
361.001	Double (McIndoe) Skin Hook	01 No
478.005-6	Needle Holder - Deep Serration 5", 6" -1 each	02 Nos
478.305-6	Needle Holder - Ryder 5", 6" -1 each	02 Nos
477.T4C-5C	Scissors - Steven's Tenotomy 4", 5" Curved - 1 each	02 Nos
477.T6C-7C	Scissors - Steven's Tenotomy 6", 7" Curved - 1 each	02 Nos
481.05P	Dissecting Forceps - Adson 5" Plain	02 Nos
481.05T	Dissecting Forceps - Adson 5" Tooth	02 Nos

Ref.	Instrument Contents	Qty
294.001	Kilner Skin Retractor Double Ended	02 Nos
480.04C-5C	Mosquito Artery Forceps 4", 5" Curved - 1 each	02 Nos
479.06C	Crile Artery Forceps - 6" Curved	02 Nos
477.S01	Scissors - Northbent Suture Cutting	01 No
477.H01	Scissors-"Heath" Stitch Cutting-4"	01 No
482.004-06	Allis Tissue Holding Forceps 4", 6" -1 each	02 Nos
492.001	Container for Basic Plastic Surgery Set	01 No



Ref.	Instrument Contents	Qty
034.001	Arch Bar	4 Nos
1551.026	Wire Reel - S.S. Premium 26 G	4 Nos
482.006	Allis Tissue Holding Forceps 6"	1 No
492.003	Alu. Anodised Box - Big	1 No
298.001	Austin Tissue Retractor - Regular	1 No
486.003	B. P. Handle No. 3	2 Nos
430.001	Bone File/Rasp Double Ended	1 No
350.004	Coupland Elevator - 4mm	1 No
479.06C	Crile Artery Forceps - 6" Curved	1 No
479.06S	Crile Artery Forceps - 6" Straight	1 No
481.05P	Dissecting Forceps - Adson 5" Plain	1 No
481.05T	Dissecting Forceps - Adson 5" Tooth	1 No
394.001	Ferguson Mouth Gag	1 No
443.001	Fickling Dental Mallet	1 No
374.004	Fine Chisels St. 4mm	1 No
341.001	Howarth Elevator	2 Nos
294.001	Kilner Skin Retractor - D/e	1 No
386.001	Lucas Currettes	1 No
351.001	Molt Elevator	2 Nos

Ref.	Instrument Contents	Qty
480.05C	Mosquito Artery Forceps 5" Curved	1 No
480.05S	Mosquito Artery Forceps 5" Straight	1 No
460.001	Mouth Props Set Of 3	1 No
478.306	Needle Holders Ryder 6"	2 Nos
324.003	Plastic Cheek Retractor - Big	1 No
324.001	Plastic Cheek Retractor - Small	1 No
273.006	Reduction Bone Holding Forceps Pointed 6"	1 No
477.I4C	Scissors - Irish Type 4.5" Curved	1 No
477.I4S	Scissors - Irish Type 4.5" Straight	1 No
477.Y6C	Scissors - Mayo 6.5" Curved	1 No
477.Y6S	Scissors - Mayo 6.5" Straight	1 No
477.M6C	Scissors - Metzenbum 6" Curved	1 No
477.M6S	Scissors - Metzenbum 6" Straight	1 No
477.S01	Scissors - Northbent Suture Cutting	2 Nos
457.002	Tongue Depressor - Serrated	1 No
485.005	Towel Forceps (Guarded) 5"	4 Nos
337.001	Warwick James Elevator - Straight	1 No
450.001	Wire Cutter - Regular	1 No
446.001	Wire Twisting Forceps (Thick)	2 Nos

STORAGE

Surgical instruments should be stored in sterilizing trays of proper size lined with soft silicone mats. Instruments should not touch with each other. We recommend using protective tips made of soft silicone tubing of the proper size and thickness. Do not use rubber or plastic protective tips, which can melt during autoclaving and cause damage to instruments.

INTENDED USE

The instruments shall be used for its intended purpose only failing of which can lead to permanent damage or breakage of the instrument.

INSPECTION

Handle the instruments with care and be sure to inspect every instrument at the end of your surgical day. Please conduct this inspection under a microscope or with magnification lens for fine instrument.

CLEANING

Manual Cleaning – All the new instruments must be decontaminated, lubricated and sterilized prior to each use. Use a mild soap solution approved for surgical instruments with a soft brush or instrument wipe to clean instrument and remove stains. Never use abrasive powders or steel wool to remove stubborn stains—these can damage the super fine finish of an instrument and can actually help cause corrosion of stainless instruments. Use distilled and treated water to rinse instruments thoroughly then carefully dry them with a hot air blower or lint-free cloth.

Ultrasonic Cleaning – This is the best method of cleaning surgical instruments. Use only cleaning solutions approved for surgical instruments. Change the solution frequently, always use treated or distilled water. Heat water to 150° Fahrenheit (65 ° C) using an immersion heater, or in a separate stainless container, if your unit does not have an automatic built-in water heater. Be sure that you and your staff members are completely familiar with the manufacturer's instructions, which came with your ultrasonic cleaner. Instruments should never touch each other in the ultrasonic cleaner. After removing instruments from the ultrasonic cleaner, rinse the instruments under running water first, and then rinse again in clean demineralized water. Instruments should never be air dried in room temperature because this will cause corrosion on stainless instruments especially in knurled and laser engraved areas. This is why it is highly recommend using a hot air blower to thoroughly dry the instruments. Do not use acidic or alkaline solutions for cleaning because it removes the instruments' protective barrier of chromium oxide and can lead to corrosion, pitting or even breakage also.

CAUTION

Keep the instruments out of touch of following solutions: Aqua Regia- Lodine – Ferric Chloride – Sulfuric Acid – Hydrochloride Acid. Any kind of corrosion will lead to rust on steel and rust particles can be transferred from one instruments to other. Therefore any corroding instrument should be removed from service to prevent the formation of rust to other instruments

LUBRICATION

It is recommended to lubricate all moving parts after each cleaning. Use only silicone or Teflon lubricants which can hold up under autoclaving. Never use ordinary lubricants or cover up even the slightest corrosion with a lubricant.

STERILIZATION

Stainless Steel instruments can be sterilized via steam autoclaving, chemical disinfectants, ethylene oxide gas, or even dry hot air. Gas and dry chemical sterilization are the best methods for stainless steel instruments, but they take a lengthy time period to accomplish the desired result. The most practical method of sterilization is heat or steam, which require less time. However, these methods can be damaging to delicate stainless steel instruments and instruments with Fiber or Teflon handles. Take care when sterilizing your instruments. And, once again, be sure that you and your staff have read and fully understood the instructions supplied by the manufacturer of your particular sterilizer. It is advisable to sterilize the instruments in the open position or disassembled because steam will only sterilize the surface where it can touch directly.

Kindly follow above instruction to enjoy long service life of the instruments.